

Protocol: IMAGENB 1 6 8

LOCALISER	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	2D
	Coil Configuration	8HRBRAIN	Pulse Sequence	Gradient Echo
	Plane	3-PLANE	Imaging Options	Seq, Fast
	Series Description	LOCALISER	Acceleration Factor	1.00
	<b>SCAN RANGE</b>		<b>IMAGE ENHANCE</b>	
	FOV	24.0	Filter Choice	None
	Slice Thickness	5.0		
	Slice Spacing	0.5		
	<b>ACQ TIMING</b>		<b>GATING/TRIGGER</b>	
	Freq	256	Auto Trigger Type	Off
	Phase	128		
	Freq DIR	Unswap		
	NEX	1.00		
	Phase FOV	1.00		
	Auto Shim	Auto		
	Phase Correction	No		
	<b>FMRI</b>		<b>MULTI-PHASE</b>	
	Initial State	Control	# of Acquisition	0
	PSD Trigger	Internal	Seperate Series	0
	Slice Order	Interleaved	Mask Phase	0
	View Order	Bottom/Up	Mask Pause	0
	# of Repetitions REST	3		
	# of Repetitions ACTIVE	3		
	<b>SAT</b>		<b>ASSET</b>	
	Tag Type	None	Slice Acceleration Factor	1.00
			Phase Acceleration Factor	1.00
	<b>TRICKS</b>		<b>CONTRAST</b>	
	Pause On/Off	On	Contrast Yes/No	No
	Auto Subtract	0		
	Auto SCIC	Off		
	<b>MULTI-STATION</b>			
	Number of Stations	0		
	Station Number	0		

Protocol: IMAGENB 1 6 8

ASSET Cal	PATIENT POSITION		IMAGING PARAMETERS		ASSET Cal
	Patient Entry	Head First	Gradient Mode	WHOLE	
	Patient Position	Supine	Imaging Mode	2D	
	Coil Configuration	8HRBRAIN	Pulse Sequence	Gradient Echo	
	Plane	AXIAL	Imaging Options	Fast, Calib	
	Series Description	ASSET Cal			
	SCAN TIMING		SCAN RANGE		
	Number of Echoes	1	FOV	30.0	
			Slice Thickness	8.0	
			Slice Spacing	0.0	
	IMAGE ENHANCE		ACQ TIMING		
	Filter Choice	None	Freq DIR	R/L	
			Auto Shim	Auto	
			Phase Correction	No	
	GATING/TRIGGER		MULTI-PHASE		
	Auto Trigger Type	Off	Seperate Series	0	
			Mask Phase	0	
			Mask Pause	0	
	SAT		TRICKS		
	Tag Type	None	Pause On/Off	On	
		Auto Subtract	0		
		Auto SCIC	Off		
CONTRAST		MULTI-STATION			
Contrast Yes/No	No	Number of Stations	0		
		Station Number	0		

ADNI MPRAGE	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	3D
	Coil Configuration	8HRBRAIN	Pulse Sequence	SPGR
	Plane	SAGITTAL	Imaging Options	EDR, Fast, IrP
	Series Description	ADNI MPRAGE	PSD Name	/usr/g/bin/adni14m4
	<b>SCAN TIMING</b>		Acceleration Factor	1.00
	<b>SCAN RANGE</b>		FOV	28.0
	Flip Angle	8	Slice Thickness	1.1
	TE	Min Full	Location per Slab	170
	Number of Echoes	1	Overlap Locations	0
	TI	900	<b>ACQ TIMING</b>	
	Receiver Bandwidth	31.25	Freq	256
	<b>IMAGE ENHANCE</b>		Phase	256
	Filter Choice	None	Freq DIR	S/I
	<b>GATING/TRIGGER</b>		NEX	1.00
	Auto Trigger Type	Off	Phase FOV	0.94
	<b>FMRI</b>		Auto Shim	Auto
	PSD Trigger	Internal	Phase Correction	No
	Slice Order	Interleaved	<b>USER CVS</b>	
	View Order	Bottom/Up	User CV6	1.00
	# of Repetitions REST	0	User CV17	1.00
	# of Repetitions ACTIVE	0	User CV18	1.00
	<b>SAT</b>		User CV23	100.00
	Tag Type	None	User CV24	2300.00
	<b>TRICKS</b>		<b>MULTI-PHASE</b>	
	Pause On/Off	On	# of Acquisition	0
	Auto Subtract	0	Seperate Series	0
	Auto SCIC	Off	Mask Phase	0
	<b>MULTI-STATION</b>		Mask Pause	0
	Number of Stations	0	<b>ASSET</b>	
			Slice Acceleration Factor	1.00
			Phase Acceleration Factor	1.00
			<b>CONTRAST</b>	
			Contrast Yes/No	No

---

Protocol: IMAGENB 1 6 8

<b>Station Number</b>
-----------------------

0

## Protocol: IMAGENB 1 6 8

T2

PATIENT POSITION	
Patient Entry	Head First
Patient Position	Supine
Coil Configuration	8HRBRAIN
Plane	OBLIQUE
Series Description	T2
SCAN TIMING	
TE	60.0
Number of Echoes	1
TR	4380.0
Echo Train Length	19
Receiver Bandwidth	25.00
IMAGE ENHANCE	
Filter Choice	None

Protocol: IMAGENB 1 6 8

FLAIR	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	2D
	Coil Configuration	8HRBRAIN	Pulse Sequence	T2flair
	Plane	OBLIQUE	Imaging Options	FC, TRF, Fast, ZIP512
	Series Description	FLAIR	Acceleration Factor	1.00
	<b>SCAN TIMING</b>		<b>SCAN RANGE</b>	
	TE	120.0	FOV	24.0
	Number of Echoes	1	Slice Thickness	4.0
	TR	8000.0	Slice Spacing	0.0
	TI	2000		
	Receiver Bandwidth	25.00		
	<b>IMAGE ENHANCE</b>		<b>ACQ TIMING</b>	
	Filter Choice	None	Freq	256
			Phase	128
			Freq DIR	A/P
			NEX	1.00
			Auto Shim	Auto
			Phase Correction	No
			Flow Direction	Freq
			Compensation	
	<b>GATING/TRIGGER</b>		<b>USER CVS</b>	
	Auto Trigger Type	Off	User CV3	2.00
			User CV20	0.0
	<b>FMRI</b>		<b>MULTI-PHASE</b>	
	PSD Trigger	Internal	# of Acquisition	0
	Slice Order	Interleaved	Seperate Series	0
	View Order	Bottom/Up	Mask Phase	0
	# of Repetitions REST	0	Mask Pause	0
	# of Repetitions ACTIVE	0		
	<b>SAT</b>		<b>ASSET</b>	
	SAT Location	I	Slice Acceleration Factor	1.00
	Tag Type	None	Phase Acceleration Factor	1.00
	<b>TRICKS</b>		<b>CONTRAST</b>	
	Pause On/Off	On	Contrast Yes/No	No
	Auto Subtract	0		
	Auto SCIC	Off		
	<b>MULTI-STATION</b>			
	Number of Stations	0		
	Station Number	0		

EPI MID (IMAGEN 2)	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	2D
	Coil Configuration	8HRBRAIN	Pulse Sequence	Gradient Echo
	Plane	OBLIQUE	Imaging Options	EPI, Asset, FMRI
	Series Description	EPI MID (IMAGEN 2)	Acceleration Factor	1.00
	<b>SCAN TIMING</b>		<b>SCAN RANGE</b>	
	Flip Angle	75	FOV	21.8
	TE	30.0	Slice Thickness	2.4
	Number of Echoes	1	Slice Spacing	1.0
	TR	2200.0		
	Number of Shots	1		
	<b>IMAGE ENHANCE</b>		<b>ACQ TIMING</b>	
	Filter Choice	None	Freq	64
			Phase	64
			Freq DIR	R/L
			NEX	1.00
			Phase FOV	1.00
			Auto Shim	Auto
			Phase Correction	Yes
	<b>GATING/TRIGGER</b>		<b>USER CVS</b>	
	Auto Trigger Type	Off	User CV0	1.00
	<b>FMRI</b>		<b>MULTI-PHASE</b>	
	Brain Wave Real Time	1	Slice per Location	191
	Paradigm String	IMAGENB - MID	Delay after Acquisition	0
	Paradigm UID	1.2.840.113819.3.6034287119 6.1192025176.127	# of Acquisition	0
	Initial State	Control	Seperate Series	0
	PSD Trigger	Internal	Delay after Acquisition without AV	0
	Slice Order	Sequential	Mask Phase	0
	View Order	Bottom/Up	Mask Pause	0
	# of Repetitions REST	1		
	# of Repetitions ACTIVE	1		
	# of Dummy Acquisition	3		
	<b>SAT</b>		<b>ASSET</b>	
	Tag Type	None	Slice Acceleration Factor	1.00
			Phase Acceleration Factor	2.00
	<b>TRICKS</b>		<b>CONTRAST</b>	
	Pause On/Off	On	Contrast Yes/No	No
	Auto Subtract	0		

---

Protocol: IMAGENB 1 6 8

Auto SCIC	Off
MULTI-STATION	
Number of Stations	0
Station Number	0



EPI stop signal (IMAGEN 2)	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	2D
	Coil Configuration	8HRBRAIN	Pulse Sequence	Gradient Echo
	Plane	OBLIQUE	Imaging Options	EPI, Asset, FMRI
	Series Description	EPI stop signal (IMAGEN 2)	Acceleration Factor	1.00
	<b>SCAN TIMING</b>		<b>SCAN RANGE</b>	
	Flip Angle	75	FOV	21.8
	TE	30.0	Slice Thickness	2.4
	TR	2200.0	Slice Spacing	1.0
	Number of Shots	1	<b>ACQ TIMING</b>	
	<b>IMAGE ENHANCE</b>		Freq	64
	Filter Choice	None	Phase	64
	<b>GATING/TRIGGER</b>		Freq DIR	R/L
	Auto Trigger Type	Off	NEX	1.00
	<b>FMRI</b>		Phase FOV	1.00
	Brain Wave Real Time	1	Auto Shim	Auto
	Paradigm String	IMAGENB - Stop Signal	Phase Correction	Yes
	Paradigm UID	1.2.840.113819.3.6034287119 6.1208955039.971	<b>USER CVS</b>	
	Initial State	Control	User CV0	1.00
	PSD Trigger	Internal	<b>MULTI-PHASE</b>	
	Slice Order	Sequential	Slice per Location	350
	View Order	Bottom/Up	Delay after Acquisition	0
	# of Repetitions REST	1	# of Acquisition	0
	# of Repetitions ACTIVE	1	Seperate Series	0
	# of Dummy Acquisition	3	Delay after Acquisition without AV	0
	<b>SAT</b>		Mask Phase	0
	Tag Type	None	Mask Pause	0
	<b>TRICKS</b>		<b>ASSET</b>	
	Pause On/Off	On	Slice Acceleration Factor	1.00
	Auto Subtract	0	Phase Acceleration Factor	2.00
	Auto SCIC	Off	<b>CONTRAST</b>	
	<b>MULTI-STATION</b>		Contrast Yes/No	No

---

Protocol: IMAGENB 1 6 8

<b>Number of Stations</b>	0
<b>Station Number</b>	0

EPI faces (IMAGEN 2)	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	2D
	Coil Configuration	8HRBRAIN	Pulse Sequence	Gradient Echo
	Plane	OBLIQUE	Imaging Options	EPI, Asset, FMRI
	Series Description	EPI faces (IMAGEN 2)	Acceleration Factor	1.00
	SCAN TIMING		SCAN RANGE	
	Flip Angle	75	FOV	21.8
	TE	30.0	Slice Thickness	2.4
	TR	2200.0	Slice Spacing	1.0
	Number of Shots	1	ACQ TIMING	
	IMAGE ENHANCE		Freq	64
	Filter Choice	None	Phase	64
	GATING/TRIGGER		Freq DIR	R/L
	Auto Trigger Type	Off	NEX	1.00
	FMRI		Phase FOV	1.00
	Brain Wave Real Time	1	Auto Shim	Auto
	Paradigm String	IMAGENB - Faces	Phase Correction	Yes
	Paradigm UID	1.2.840.113819.3.6034287119 6.1208954758.406	USER CVS	
	Initial State	Control	User CV0	1.00
	PSD Trigger	Internal	MULTI-PHASE	
	Slice Order	Sequential	Slice per Location	197
	View Order	Bottom/Up	Delay after Acquisition	0
	# of Repetitions REST	1	# of Acquisition	0
	# of Repetitions ACTIVE	1	Seperate Series	0
	# of Dummy Acquisition	3	Delay after Acquisition without AV	0
	SAT		Mask Phase	0
	Tag Type	None	Mask Pause	0
	TRICKS		ASSET	
	Pause On/Off	On	Slice Acceleration Factor	1.00
	Auto Subtract	0	Phase Acceleration Factor	2.00
	Auto SCIC	Off	CONTRAST	
	MULTI-STATION		Contrast Yes/No	No

---

Protocol: IMAGENB 1 6 8

<b>Number of Stations</b>	0
<b>Station Number</b>	0

FIELD MAP 8 CHANNEL	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	2D
	Coil Configuration	8HRBRAIN	Pulse Sequence	Spin Echo
	Plane	OBLIQUE	Imaging Options	FC
	Series Description	FIELD MAP 8 CHANNEL	PSD Name	/usr/g/cns/spgtdjl/fcmempase
	SCAN TIMING		Acceleration Factor	1.00
	TE	16.0	SCAN RANGE	
	Number of Echoes	1	FOV	25.6
	TR	2200.0	Slice Thickness	4.0
	Receiver Bandwidth	8.06	Slice Spacing	0.0
	IMAGE ENHANCE		ACQ TIMING	
	Filter Choice	None	Freq	64
	GATING/TRIGGER		Phase	64
	Auto Trigger Type	Off	Freq DIR	A/P
	FMRI		NEX	1.00
	Brain Wave Real Time	1	Phase FOV	1.00
	Paradigm String	IMAGEN - Global_correc	Auto Shim	Auto
	Paradigm UID	1.2.840.113819.3.6034287119 6.1192025176.127	Phase Correction	No
	Initial State	Control	USER CVS	
	PSD Trigger	Internal	User CV0	1136.00
	Slice Order	Sequential	MULTI-PHASE	
	View Order	Bottom/Up	# of Acquisition	0
	# of Repetitions REST	1	Seperate Series	0
	# of Repetitions ACTIVE	1	Mask Phase	0
	# of Dummy Acquisition	3	Mask Pause	0
	SAT		ASSET	
	Tag Type	None	Slice Acceleration Factor	1.00
	TRICKS		Phase Acceleration Factor	1.00
	Pause On/Off	On	CONTRAST	
	Auto Subtract	0	Contrast Yes/No	No
	Auto SCIC	Off	MULTI-STATION	

---

Protocol: IMAGENB 1 6 8

<b>Number of Stations</b>	0
<b>Station Number</b>	0

EPI REST (IMAGEN 2)	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	2D
	Coil Configuration	8HRBRAIN	Pulse Sequence	Gradient Echo
	Plane	OBLIQUE	Imaging Options	EPI, Asset, FMRI
	Series Description	EPI REST (IMAGEN 2)	Acceleration Factor	1.00
	<b>SCAN TIMING</b>		<b>SCAN RANGE</b>	
	Flip Angle	75	FOV	21.8
	TE	30.0	Slice Thickness	2.4
	Number of Echoes	1	Slice Spacing	1.0
	TR	2200.0		
	Number of Shots	1		
	<b>IMAGE ENHANCE</b>		<b>ACQ TIMING</b>	
	Filter Choice	None	Freq	64
			Phase	64
			Freq DIR	R/L
			NEX	1.00
			Phase FOV	1.00
			Auto Shim	Auto
			Phase Correction	Yes
	<b>GATING/TRIGGER</b>		<b>USER CVS</b>	
	Auto Trigger Type	Off	User CV0	1.00
	<b>FMRI</b>		<b>MULTI-PHASE</b>	
EPI REST (IMAGEN 2)	Brain Wave Real Time	1	Slice per Location	164
	Paradigm String	IMAGENB - RESTING STATE	Delay after Acquisition	0
	Paradigm UID	1.2.840.113819.3.6034287119 6.1224159165.89	# of Acquisition	0
	Initial State	Control	Seperate Series	0
	PSD Trigger	Internal	Delay after Acquisition without AV	0
	Slice Order	Sequential	Mask Phase	0
	View Order	Bottom/Up	Mask Pause	0
	# of Repetitions REST	1		
	# of Repetitions ACTIVE	1		
	# of Dummy Acquisition	3		
	<b>SAT</b>		<b>ASSET</b>	
	Tag Type	None	Slice Acceleration Factor	1.00
			Phase Acceleration Factor	2.00
	<b>TRICKS</b>		<b>CONTRAST</b>	
	Pause On/Off	On	Contrast Yes/No	No
	Auto Subtract	0		
	Auto SCIC	Off		

---

Protocol: IMAGENB 1 6 8

**MULTI-STATION**

**Number of Stations**      0

**Station Number**        0



PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Gradient Mode	ZOOM
Patient Position	Supine	Imaging Mode	2D
Coil Configuration	8HRBRAIN	Pulse Sequence	Spin Echo
Plane	OBLIQUE	Imaging Options	Gat, EPI, DIFF, Asset
Series Description	DTI	PSD Name	/usr/g/cns/bin/dti_CNS_08
		Acceleration Factor	1.00
SCAN TIMING		SCAN RANGE	
TE	Minimum	FOV	30.7
Number of Echoes	1	Slice Thickness	2.4
Number of Shots	1	Slice Spacing	0.0
IMAGE ENHANCE		ACQ TIMING	
Filter Choice	None	Freq	128
		Phase	128
		Freq DIR	R/L
		NEX	1.00
		Phase FOV	1.00
		Auto Shim	Auto
		Phase Correction	Yes
GATING/TRIGGER		USER CVS	
Auto Trigger Type	Off	User CV11	1.00
Trigger Type	5		
Minimum Trigger Delay	1		
Phase Select	Single Phase		
Cardiac Phase	1		
Cardiac Slices	1		
FMRI		MULTI-PHASE	
PSD Trigger	Internal	# of Acquisition	0
Slice Order	Interleaved	Seperate Series	0
View Order	Bottom/Up	Inter-Delay	Minimum
# of Repetitions REST	0	Mask Phase	0
# of Repetitions ACTIVE	0	Mask Pause	0
SAT		DIFFUSION	
Tag Type	None	B Value	1300
		Optimized TE	No
		Diffusion Directions	Tensor
		Number of Diffusion Directions	32
		Number of T2 Images	4
		Dual Spin Echo	On
		Diffusion Tenser	No Selection
		Processing Output	

---

Protocol: IMAGENB 1 6 8

<b>ASSET</b>		Recon All Images	<i>Off</i>
Slice Acceleration Factor		1.00	
Phase Acceleration Factor		2.00	
<b>CONTRAST</b>		<b>TRICKS</b>	
Contrast Yes/No		No	
		# of RR Intervals	30
		Trigger Delay	15
		Pause On/Off	On
		Auto Subtract	0
		Auto SCIC	Off
		<b>MULTI-STATION</b>	
		Number of Stations	0
		Station Number	0