

Protocol: IMAGEN_SESSION_1

LOCALISER	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	2D
	Coil Configuration	8HRBRAIN	Pulse Sequence	Gradient Echo
	Plane	3-PLANE	Imaging Options	Seq, Fast
	Series Description	LOCALISER	Acceleration Factor	1.00
	SCAN TIMING		SCAN RANGE	
	Number of Echoes	1	FOV	24.0
			Slice Thickness	5.0
			Slice Spacing	0.5
	IMAGE ENHANCE		ACQ TIMING	
	Filter Choice	None	Freq	256
			Phase	128
			Freq DIR	Unswap
			NEX	1.00
			Phase FOV	1.00
			Auto Shim	Auto
			Phase Correction	No
	GATING/TRIGGER		FMRI	
	Auto Trigger Type	Off	Initial State	Control
			PSD Trigger	Internal
			Slice Order	Interleaved
			View Order	Bottom/Up
			# of Repetitions REST	3
			# of Repetitions ACTIVE	3
	MULTI-PHASE		SAT	
	# of Acquisition	0	Tag Type	None
	Seperate Series	0		
	Mask Phase	0		
	Mask Pause	0		
	ASSET		TRICKS	
	Slice Acceleration Factor	1.00	Pause On/Off	On
	Phase Acceleration Factor	1.00	Auto Subtract	0
			Auto SCIC	Off
	CONTRAST		MULTI-STATION	
	Contrast Yes/No	No	Number of Stations	0
			Station Number	0

Protocol: IMAGEN_SESSION_1

ASSET Cal	PATIENT POSITION		IMAGING PARAMETERS		ASSET Cal
	Patient Entry	Head First	Gradient Mode	WHOLE	
	Patient Position	Supine	Imaging Mode	2D	
	Coil Configuration	8HRBRAIN	Pulse Sequence	Gradient Echo	
	Plane	AXIAL	Imaging Options	Fast, Calib	
	Series Description	ASSET Cal			
	SCAN TIMING		SCAN RANGE		
	Number of Echoes	1	FOV	30.0	
			Slice Thickness	8.0	
			Slice Spacing	0.0	
	IMAGE ENHANCE		ACQ TIMING		
	Filter Choice	None	Freq DIR	R/L	
			Auto Shim	Auto	
			Phase Correction	No	
	GATING/TRIGGER		MULTI-PHASE		
	Auto Trigger Type	Off	Seperate Series	0	
			Mask Phase	0	
			Mask Pause	0	
SAT		TRICKS			
Tag Type	None	Pause On/Off	On		
		Auto Subtract	0		
		Auto SCIC	Off		
CONTRAST		MULTI-STATION			
Contrast Yes/No	No	Number of Stations	0		
		Station Number	0		

Protocol: IMAGEN_SESSION_1

T2

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Gradient Mode	ZOOM
Patient Position	Supine	Imaging Mode	2D
Coil Configuration	8HRBRAIN	Pulse Sequence	FRFSE-XL
Plane	OBLIQUE	Imaging Options	FC, TRF, Fast, ZIP512, Asset, FR
Series Description	T2	Acceleration Factor	1.00
SCAN TIMING		SCAN RANGE	
TE	60.0	FOV	24.0
Number of Echoes	1	Slice Thickness	4.0
TR	4380.0	Slice Spacing	0.0
Echo Train Length	19		
Receiver Bandwidth	25.00		
IMAGE ENHANCE		ACQ TIMING	
Filter Choice	None	Freq	320
		Phase	320
		Freq DIR	A/P
		NEX	1.00
		# of Acq. Before Pause	0
		Phase FOV	0.75
		Auto Shim	Auto
		Phase Correction	No
		Flow Direction	Freq
		Compensation	
GATING/TRIGGER		USER CVS	
Auto Trigger Type	Off	User CV20	0.0
FMRI		MULTI-PHASE	
PSD Trigger	Internal	# of Acquisition	0
Slice Order	Interleaved	Seperate Series	0
View Order	Bottom/Up	Mask Phase	0
# of Repetitions REST	0	Mask Pause	0
# of Repetitions ACTIVE	0		
SAT		ASSET	
SAT Location	I	Slice Acceleration Factor	1.00
Tag Type	None	Phase Acceleration Factor	2.00
TRICKS		CONTRAST	
Pause On/Off	On	Contrast Yes/No	No
Auto Subtract	0		
Auto SCIC	Off		
MULTI-STATION			
Number of Stations	0		
Station Number	0		

T2

Protocol: IMAGEN_SESSION_1

FLAIR	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	2D
	Coil Configuration	8HRBRAIN	Pulse Sequence	T2flair
	Plane	OBLIQUE	Imaging Options	FC, TRF, Fast, ZIP512
	Series Description	FLAIR	Acceleration Factor	1.00
	SCAN TIMING		SCAN RANGE	
	TE	120.0	FOV	24.0
	Number of Echoes	1	Slice Thickness	4.0
	TR	8000.0	Slice Spacing	0.0
	TI	2000		
	Receiver Bandwidth	25.00		
	IMAGE ENHANCE		ACQ TIMING	
	Filter Choice	None	Freq	256
			Phase	128
			Freq DIR	A/P
			NEX	1.00
			Auto Shim	Auto
			Phase Correction	No
			Flow Direction	Freq
			Compensation	
	GATING/TRIGGER		USER CVS	
	Auto Trigger Type	Off	User CV3	2.00
			User CV20	0.0
	FMRI		MULTI-PHASE	
	PSD Trigger	Internal	# of Acquisition	0
	Slice Order	Interleaved	Seperate Series	0
	View Order	Bottom/Up	Mask Phase	0
	# of Repetitions REST	0	Mask Pause	0
	# of Repetitions ACTIVE	0		
	SAT		ASSET	
	SAT Location	I	Slice Acceleration Factor	1.00
	Tag Type	None	Phase Acceleration Factor	1.00
	TRICKS		CONTRAST	
	Pause On/Off	On	Contrast Yes/No	No
	Auto Subtract	0		
	Auto SCIC	Off		
	MULTI-STATION			
	Number of Stations	0		
	Station Number	0		

EPI faces	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	2D
	Coil Configuration	8HRBRAIN	Pulse Sequence	Gradient Echo
	Plane	OBLIQUE	Imaging Options	EPI, Asset, FMRI
	Series Description	EPI faces	Acceleration Factor	1.00
	SCAN TIMING		SCAN RANGE	
	Flip Angle	75	FOV	21.8
	TE	30.0	Slice Thickness	2.4
	Number of Echoes	1	Slice Spacing	1.0
	TR	2200.0		
	Number of Shots	1		
	IMAGE ENHANCE		ACQ TIMING	
	Filter Choice	None	Freq	64
			Phase	64
			Freq DIR	R/L
			NEX	1.00
			Phase FOV	1.00
			Auto Shim	Auto
			Phase Correction	Yes
	GATING/TRIGGER		USER CVS	
	Auto Trigger Type	Off	User CV0	1.00
	FMRI		MULTI-PHASE	
	Brain Wave Real Time	1	Slice per Location	160
	Paradigm String	IMAGEN - Faces	Delay after Acquisition	0
	Paradigm UID	1.2.840.113819.3.6034287119 6.1208954758.406	# of Acquisition	0
	Initial State	Control	Seperate Series	0
	PSD Trigger	Internal	Delay after Acquisition without AV	0
	Slice Order	Sequential	Mask Phase	0
	View Order	Bottom/Up	Mask Pause	0
	# of Repetitions REST	1		
	# of Repetitions ACTIVE	1		
	# of Dummy Acquisition	3		
	SAT		ASSET	
	Tag Type	None	Slice Acceleration Factor	1.00
			Phase Acceleration Factor	2.00
	TRICKS		CONTRAST	
	Pause On/Off	On	Contrast Yes/No	No
	Auto Subtract	0		

Protocol: IMAGEN_SESSION_1

Auto SCIC	Off
MULTI-STATION	
Number of Stations	0
Station Number	0

EPI stop signal	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	2D
	Coil Configuration	8HRBRAIN	Pulse Sequence	Gradient Echo
	Plane	OBLIQUE	Imaging Options	EPI, Asset, FMRI
	Series Description	EPI stop signal	Acceleration Factor	1.00
	SCAN TIMING		SCAN RANGE	
	Flip Angle	75	FOV	21.8
	TE	30.0	Slice Thickness	2.4
	Number of Echoes	1	Slice Spacing	1.0
	TR	2200.0		
	Number of Shots	1		
	IMAGE ENHANCE		ACQ TIMING	
	Filter Choice	None	Freq	64
			Phase	64
			Freq DIR	R/L
			NEX	1.00
			Phase FOV	1.00
			Auto Shim	Auto
			Phase Correction	Yes
	GATING/TRIGGER		USER CVS	
	Auto Trigger Type	Off	User CV0	1.00
	FMRI		MULTI-PHASE	
	Brain Wave Real Time	1	Slice per Location	444
	Paradigm String	IMAGEN - Stop Signal	Delay after Acquisition	0
	Paradigm UID	1.2.840.113819.3.6034287119 6.1208955039.971	# of Acquisition	0
	Initial State	Control	Seperate Series	0
	PSD Trigger	Internal	Delay after Acquisition without AV	0
	Slice Order	Sequential	Mask Phase	0
	View Order	Bottom/Up	Mask Pause	0
	# of Repetitions REST	1		
	# of Repetitions ACTIVE	1		
	# of Dummy Acquisition	3		
	SAT		ASSET	
	Tag Type	None	Slice Acceleration Factor	1.00
			Phase Acceleration Factor	2.00
	TRICKS		CONTRAST	
	Pause On/Off	On	Contrast Yes/No	No
	Auto Subtract	0		
	Auto SCIC	Off		

Protocol: IMAGEN_SESSION_1

MULTI-STATION

Number of Stations 0

Station Number 0

FIELD MAP 8 CHANNEL	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	2D
	Coil Configuration	8HRBRAIN	Pulse Sequence	Spin Echo
	Plane	OBLIQUE	Imaging Options	FC
	Series Description	FIELD MAP 8 CHANNEL	PSD Name	/usr/g/cns/spgtdjl/fcmempase
	SCAN TIMING		Acceleration Factor	1.00
	TE	16.0	SCAN RANGE	
	Number of Echoes	1	FOV	25.6
	TR	2200.0	Slice Thickness	4.0
	Receiver Bandwidth	8.06	Slice Spacing	0.0
	IMAGE ENHANCE		ACQ TIMING	
	Filter Choice	None	Freq	64
	GATING/TRIGGER		Phase	64
	Auto Trigger Type	Off	Freq DIR	A/P
	FMRI		NEX	1.00
	Brain Wave Real Time	1	Phase FOV	1.00
	Paradigm String	IMAGEN - Global_correc	Auto Shim	Auto
	Paradigm UID	1.2.840.113819.3.6034287119 6.1192025176.127	Phase Correction	No
	Initial State	Control	USER CVS	
	PSD Trigger	Internal	User CV0	1136.00
	Slice Order	Sequential	MULTI-PHASE	
	View Order	Bottom/Up	# of Acquisition	0
	# of Repetitions REST	1	Seperate Series	0
	# of Repetitions ACTIVE	1	Mask Phase	0
	# of Dummy Acquisition	3	Mask Pause	0
	SAT		ASSET	
	Tag Type	None	Slice Acceleration Factor	1.00
	TRICKS		Phase Acceleration Factor	1.00
	Pause On/Off	On	CONTRAST	
	Auto Subtract	0	Contrast Yes/No	No
	Auto SCIC	Off	MULTI-STATION	

Protocol: IMAGEN_SESSION_1

Number of Stations	0
Station Number	0

ADNI MPRAGE	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	3D
	Coil Configuration	8HRBRAIN	Pulse Sequence	SPGR
	Plane	SAGITTAL	Imaging Options	EDR, Fast, IrP
	Series Description	ADNI MPRAGE	PSD Name	/usr/g/bin/adni14m4
	SCAN TIMING		Acceleration Factor	1.00
	Flip Angle	8	SCAN RANGE	
	TE	Min Full	FOV	28.0
	Number of Echoes	1	Slice Thickness	1.1
	TI	900	Location per Slab	170
	Receiver Bandwidth	31.25	Overlap Locations	0
	IMAGE ENHANCE		ACQ TIMING	
	Filter Choice	None	Freq	256
			Phase	256
			Freq DIR	S/I
			NEX	1.00
			Phase FOV	0.94
			Auto Shim	Auto
	GATING/TRIGGER		Phase Correction	No
	Auto Trigger Type	Off	USER CVS	
			User CV6	1.00
			User CV17	1.00
			User CV18	1.00
			User CV23	100.00
			User CV24	2300.00
	FMRI		MULTI-PHASE	
	PSD Trigger	Internal	# of Acquisition	0
	Slice Order	Interleaved	Seperate Series	0
	View Order	Bottom/Up	Mask Phase	0
	# of Repetitions REST	0	Mask Pause	0
	# of Repetitions ACTIVE	0	ASSET	
	SAT		Slice Acceleration Factor	1.00
	Tag Type	None	Phase Acceleration Factor	1.00
	TRICKS		CONTRAST	
	Pause On/Off	On	Contrast Yes/No	No
	Auto Subtract	0		
	Auto SCIC	Off		
	MULTI-STATION			
	Number of Stations	0		

Protocol: IMAGEN_SESSION_1

Station Number

0