

Protocol: IMAGEN\_SESSION\_2

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Gradient Mode	ZOOM
Patient Position	Supine	Imaging Mode	2D
Coil Configuration	8HRBRAIN	Pulse Sequence	Gradient Echo
Plane	3-PLANE	Imaging Options	Seq, Fast
Series Description	LOCALISER	Acceleration Factor	1.00
SCAN TIMING		SCAN RANGE	
Number of Echoes	1	FOV	24.0
		Slice Thickness	5.0
		Slice Spacing	0.5
IMAGE ENHANCE		ACQ TIMING	
Filter Choice	None	Freq	256
		Phase	128
		Freq DIR	Unswap
		NEX	1.00
		Phase FOV	1.00
		Auto Shim	Auto
		Phase Correction	No
GATING/TRIGGER		FMRI	
Auto Trigger Type	Off	PSD Trigger	Internal
		Slice Order	Interleaved
		View Order	Bottom/Up
		# of Repetitions REST	0
		# of Repetitions ACTIVE	0
MULTI-PHASE		SAT	
# of Acquisition	0	Tag Type	None
Seperate Series	0		
Mask Phase	0		
Mask Pause	0		
ASSET		TRICKS	
Slice Acceleration Factor	1.00	Pause On/Off	On
Phase Acceleration Factor	1.00	Auto Subtract	0
		Auto SCIC	Off
CONTRAST		MULTI-STATION	
Contrast Yes/No	No	Number of Stations	0
		Station Number	0

## Protocol: IMAGEN\_SESSION\_2

ASSET Cal	PATIENT POSITION		IMAGING PARAMETERS		ASSET Cal
	Patient Entry	Head First	Gradient Mode	WHOLE	
	Patient Position	Supine	Imaging Mode	2D	
	Coil Configuration	8HRBRAIN	Pulse Sequence	Gradient Echo	
	Plane	AXIAL	Imaging Options	Fast, Calib	
	Series Description	ASSET Cal	Acceleration Factor	1.00	
	SCAN TIMING		SCAN RANGE		
	Number of Echoes	1	FOV	30.0	
			Slice Thickness	8.0	
			Slice Spacing	0.0	
	IMAGE ENHANCE		ACQ TIMING		
	Filter Choice	None	Freq DIR	R/L	
			Auto Shim	Auto	
			Phase Correction	No	
	GATING/TRIGGER		FMRI		
	Auto Trigger Type	Off	PSD Trigger	Internal	
			Slice Order	Interleaved	
			View Order	Bottom/Up	
			# of Repetitions REST	0	
			# of Repetitions ACTIVE	0	
	MULTI-PHASE		SAT		
	# of Acquisition	0	Tag Type	None	
	Seperate Series	0			
	Mask Phase	0			
	Mask Pause	0			
	ASSET		TRICKS		
	Slice Acceleration Factor	1.00	Pause On/Off	On	
	Phase Acceleration Factor	1.00	Auto Subtract	0	
			Auto SCIC	Off	
	CONTRAST		MULTI-STATION		
	Contrast Yes/No	No	Number of Stations	0	
			Station Number	0	

short MPRAGE	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	3D
	Coil Configuration	8HRBRAIN	Pulse Sequence	SPGR
	Plane	SAGITTAL	Imaging Options	EDR, Fast, IrP, Asset
	Series Description	short MPRAGE	PSD Name	/usr/g/bin/adni14m4
	<b>SCAN TIMING</b>		Acceleration Factor	1.00
	<b>SCAN RANGE</b>		FOV	28.0
	Flip Angle	8	Slice Thickness	1.1
	TE	Min Full	Location per Slab	170
	Number of Echoes	1	Overlap Locations	0
	TI	900	<b>ACQ TIMING</b>	
	Receiver Bandwidth	31.25	Freq	256
	<b>IMAGE ENHANCE</b>		Phase	128
	Filter Choice	None	Freq DIR	S/I
	<b>GATING/TRIGGER</b>		NEX	1.00
	Auto Trigger Type	Off	Phase FOV	0.94
	<b>FMRI</b>		Auto Shim	Auto
	PSD Trigger	Internal	Phase Correction	No
	Slice Order	Interleaved	<b>USER CVS</b>	
	View Order	Bottom/Up	User CV6	1.00
	# of Repetitions REST	0	User CV17	1.00
	# of Repetitions ACTIVE	0	User CV18	1.00
	<b>SAT</b>		User CV23	100.00
	Tag Type	None	User CV24	2300.00
	<b>TRICKS</b>		<b>MULTI-PHASE</b>	
	Pause On/Off	On	# of Acquisition	0
	Auto Subtract	0	Seperate Series	0
	Auto SCIC	Off	Mask Phase	0
	<b>MULTI-STATION</b>		Mask Pause	0
	Number of Stations	0	<b>ASSET</b>	
			Slice Acceleration Factor	1.00
			Phase Acceleration Factor	2.00
			<b>CONTRAST</b>	
			Contrast Yes/No	No

---

Protocol: IMAGEN\_SESSION\_2

Station Number	0
----------------	---

EPI short MID	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	2D
	Coil Configuration	8HRBRAIN	Pulse Sequence	Gradient Echo
	Plane	OBLIQUE	Imaging Options	EPI, Asset, FMRI
	Series Description	EPI short MID	Acceleration Factor	1.00
	<b>SCAN TIMING</b>		<b>SCAN RANGE</b>	
	Flip Angle	75	FOV	21.8
	TE	30.0	Slice Thickness	2.4
	Number of Echoes	1	Slice Spacing	1.0
	TR	2200.0		
	Number of Shots	1		
	<b>IMAGE ENHANCE</b>		<b>ACQ TIMING</b>	
	Filter Choice	None	Freq	64
			Phase	64
			Freq DIR	R/L
			NEX	1.00
			Phase FOV	1.00
			Auto Shim	Auto
			Phase Correction	Yes
	<b>GATING/TRIGGER</b>		<b>USER CVS</b>	
	Auto Trigger Type	Off	User CV0	1.00
	<b>FMRI</b>		<b>MULTI-PHASE</b>	
	Brain Wave Real Time	1	Slice per Location	300
	Paradigm String	IMAGEN - Short MID	Delay after Acquisition	0
	Paradigm UID	1.2.840.113819.3.6034287119 6.1192025176.127	# of Acquisition	0
	Initial State	Control	Seperate Series	0
	PSD Trigger	Internal	Delay after Acquisition without AV	0
	Slice Order	Sequential	Mask Phase	0
	View Order	Bottom/Up	Mask Pause	0
	# of Repetitions REST	1		
	# of Repetitions ACTIVE	1		
	# of Dummy Acquisition	3		
	<b>SAT</b>		<b>ASSET</b>	
	Tag Type	None	Slice Acceleration Factor	1.00
			Phase Acceleration Factor	2.00
	<b>TRICKS</b>		<b>CONTRAST</b>	
	Pause On/Off	On	Contrast Yes/No	No
	Auto Subtract	0		
	Auto SCIC	Off		

---

Protocol: IMAGEN\_SESSION\_2

**MULTI-STATION**

**Number of Stations**      0

**Station Number**        0

EPI global

PATIENT POSITION	
Patient Entry	Head First
Patient Position	Supine
Coil Configuration	8HRBRAIN
Plane	OBLIQUE
Series Description	EPI global
SCAN TIMING	
Flip Angle	75
TE	30.0
Number of Echoes	1
TR	2200.0
Number of Shots	1
IMAGE ENHANCE	
Filter Choice	None
GATING/TRIGGER	
Auto Trigger Type	Off
FMRI	
Brain Wave Real Time	1
Paradigm String	IMAGEN - Global
Paradigm UID	1.2.840.113819.3.6034287119 6.1192025176.127
Initial State	Control
PSD Trigger	Internal
Slice Order	Sequential
View Order	Bottom/Up
# of Repetitions REST	1
# of Repetitions ACTIVE	1
# of Dummy Acquisition	3
SAT	
Tag Type	None
TRICKS	
Pause On/Off	On
Auto Subtract	0
Auto SCIC	Off

IMAGING PARAMETERS	
Gradient Mode	ZOOM
Imaging Mode	2D
Pulse Sequence	Gradient Echo
Imaging Options	EPI, Asset, FMRI
Acceleration Factor	1.00
SCAN RANGE	
FOV	21.8
Slice Thickness	2.4
Slice Spacing	1.0
ACQ TIMING	
Freq	64
Phase	64
Freq DIR	R/L
NEX	1.00
Phase FOV	1.00
Auto Shim	Auto
Phase Correction	Yes
USER CVS	
User CV0	1.00
MULTI-PHASE	
Slice per Location	140
Delay after Acquisition	0
# of Acquisition	0
Seperate Series	0
Delay after Acquisition without AV	0
Mask Phase	0
Mask Pause	0
ASSET	
Slice Acceleration Factor	1.00
Phase Acceleration Factor	2.00
CONTRAST	
Contrast Yes/No	No

EPI global

---

Protocol: IMAGEN\_SESSION\_2

**MULTI-STATION**

**Number of Stations**      0

**Station Number**        0



EPI BH calibration	<b>PATIENT POSITION</b>		<b>IMAGING PARAMETERS</b>	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	2D
	Coil Configuration	8HRBRAIN	Pulse Sequence	Gradient Echo
	Plane	OBLIQUE	Imaging Options	EPI, Asset, FMRI
	Series Description	EPI BH calibration	Acceleration Factor	1.00
	<b>SCAN TIMING</b>		<b>SCAN RANGE</b>	
	Flip Angle	75	FOV	21.8
	TE	30.0	Slice Thickness	2.4
	Number of Echoes	1	Slice Spacing	1.0
	TR	2200.0		
	Number of Shots	1		
	<b>IMAGE ENHANCE</b>		<b>ACQ TIMING</b>	
	Filter Choice	None	Freq	64
			Phase	64
			Freq DIR	R/L
			NEX	1.00
			Phase FOV	1.00
			Auto Shim	Auto
			Phase Correction	Yes
	<b>GATING/TRIGGER</b>		<b>USER CVS</b>	
	Auto Trigger Type	Off	User CV0	1.00
	<b>FMRI</b>		<b>MULTI-PHASE</b>	
	Brain Wave Real Time	1	Slice per Location	156
	Paradigm String	IMAGEN - BH	Delay after Acquisition	0
	Paradigm UID	1.2.840.113819.3.6034287119 6.1224159165.89	# of Acquisition	0
	Initial State	Control	Seperate Series	0
	PSD Trigger	Internal	Delay after Acquisition without AV	0
	Slice Order	Sequential	Mask Phase	0
	View Order	Bottom/Up	Mask Pause	0
	# of Repetitions REST	1		
	# of Repetitions ACTIVE	1		
	# of Dummy Acquisition	3		
	<b>SAT</b>		<b>ASSET</b>	
	Tag Type	None	Slice Acceleration Factor	1.00
			Phase Acceleration Factor	2.00
	<b>TRICKS</b>		<b>CONTRAST</b>	
	Pause On/Off	On	Contrast Yes/No	No
	Auto Subtract	0		
	Auto SCIC	Off		

---

Protocol: IMAGEN\_SESSION\_2

**MULTI-STATION**

**Number of Stations**      0

**Station Number**        0

EPI rest

PATIENT POSITION	
Patient Entry	Head First
Patient Position	Supine
Coil Configuration	8HRBRAIN
Plane	OBLIQUE
Series Description	EPI rest
SCAN TIMING	
Flip Angle	75
TE	30.0
Number of Echoes	1
TR	2200.0
Number of Shots	1
IMAGE ENHANCE	
Filter Choice	None
GATING/TRIGGER	
Auto Trigger Type	Off
FMRI	
Brain Wave Real Time	1
Paradigm String	IMAGEN - RESTING STATE
Paradigm UID	1.2.840.113819.3.6034287119 6.1224159165.89
Initial State	Control
PSD Trigger	Internal
Slice Order	Sequential
View Order	Bottom/Up
# of Repetitions REST	1
# of Repetitions ACTIVE	1
# of Dummy Acquisition	3
SAT	
Tag Type	None
TRICKS	
Pause On/Off	On
Auto Subtract	0
Auto SCIC	Off

EPI rest

IMAGING PARAMETERS	
Gradient Mode	ZOOM
Imaging Mode	2D
Pulse Sequence	Gradient Echo
Imaging Options	EPI, Asset, FMRI
Acceleration Factor	1.00
SCAN RANGE	
FOV	21.8
Slice Thickness	2.4
Slice Spacing	1.0
ACQ TIMING	
Freq	64
Phase	64
Freq DIR	R/L
NEX	1.00
Phase FOV	1.00
Auto Shim	Auto
Phase Correction	Yes
USER CVS	
User CV0	1.00
MULTI-PHASE	
Slice per Location	187
Delay after Acquisition	0
# of Acquisition	0
Seperate Series	0
Delay after Acquisition without AV	0
Mask Phase	0
Mask Pause	0
ASSET	
Slice Acceleration Factor	1.00
Phase Acceleration Factor	2.00
CONTRAST	
Contrast Yes/No	No

---

Protocol: IMAGEN\_SESSION\_2

**MULTI-STATION**

**Number of Stations**      0

**Station Number**        0

FIELDMAP 8 CHANNEL	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	Head First	Gradient Mode	ZOOM
	Patient Position	Supine	Imaging Mode	2D
	Coil Configuration	8HRBRAIN	Pulse Sequence	Spin Echo
	Plane	OBLIQUE	Imaging Options	FC
	Series Description	FIELDMAP 8 CHANNEL	PSD Name	/usr/g/cns/spgtdjl/fcmempase
	SCAN TIMING		Acceleration Factor	1.00
	TE	16.0	SCAN RANGE	
	Number of Echoes	1	FOV	25.6
	TR	2200.0	Slice Thickness	4.0
	Receiver Bandwidth	8.06	Slice Spacing	0.0
	IMAGE ENHANCE		ACQ TIMING	
	Filter Choice	None	Freq	64
	GATING/TRIGGER		Phase	64
	Auto Trigger Type	Off	Freq DIR	A/P
	FMRI		NEX	1.00
	Brain Wave Real Time	1	Phase FOV	1.00
	Paradigm String	IMAGEN - Global_correc	Auto Shim	Auto
	Paradigm UID	1.2.840.113819.3.6034287119 6.1192025176.127	Phase Correction	No
	Initial State	Control	USER CVS	
	PSD Trigger	Internal	User CV0	1136.00
	Slice Order	Sequential	MULTI-PHASE	
	View Order	Bottom/Up	# of Acquisition	0
	# of Repetitions REST	1	Seperate Series	0
	# of Repetitions ACTIVE	1	Mask Phase	0
	# of Dummy Acquisition	3	Mask Pause	0
	SAT		ASSET	
	Tag Type	None	Slice Acceleration Factor	1.00
	TRICKS		Phase Acceleration Factor	1.00
	Pause On/Off	On	CONTRAST	
	Auto Subtract	0	Contrast Yes/No	No
	Auto SCIC	Off	MULTI-STATION	

---

Protocol: IMAGEN\_SESSION\_2

<b>Number of Stations</b>	0
<b>Station Number</b>	0

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Head First	Gradient Mode	ZOOM
Patient Position	Supine	Imaging Mode	2D
Coil Configuration	8HRBRAIN	Pulse Sequence	Spin Echo
Plane	OBLIQUE	Imaging Options	Gat, EPI, DIFF, Asset
Series Description	DTI	PSD Name	/usr/g/cns/bin/dti_CNS_08
		Acceleration Factor	1.00
SCAN TIMING		SCAN RANGE	
TE	Minimum	FOV	30.7
Number of Echoes	1	Slice Thickness	2.4
Number of Shots	1	Slice Spacing	0.0
IMAGE ENHANCE		ACQ TIMING	
Filter Choice	None	Freq	128
		Phase	128
		Freq DIR	R/L
		NEX	1.00
		Phase FOV	1.00
		Auto Shim	Auto
		Phase Correction	Yes
GATING/TRIGGER		USER CVS	
Auto Trigger Type	Off	User CV11	1.00
Trigger Type	5		
Minimum Trigger Delay	1		
Phase Select	Single Phase		
Cardiac Phase	1		
Cardiac Slices	1		
FMRI		MULTI-PHASE	
PSD Trigger	Internal	# of Acquisition	0
Slice Order	Interleaved	Seperate Series	0
View Order	Bottom/Up	Inter-Delay	Minimum
# of Repetitions REST	0	Mask Phase	0
# of Repetitions ACTIVE	0	Mask Pause	0
SAT		DIFFUSION	
Tag Type	None	B Value	1300
		Optimized TE	No
		Diffusion Directions	Tensor
		Number of Diffusion Directions	32
		Number of T2 Images	4
		Dual Spin Echo	On
		Diffusion Tenser	No Selection
		Processing Output	

---

Protocol: IMAGEN\_SESSION\_2

<b>ASSET</b>		<b>Recon All Images</b>	<i>Off</i>
<b>Slice Acceleration Factor</b>		<b>TRICKS</b>	
<i>1.00</i>		<b># of RR Intervals</b>	<i>30</i>
<b>Phase Acceleration Factor</b>		<b>Trigger Delay</b>	<i>15</i>
<i>2.00</i>		<b>Pause On/Off</b>	<i>On</i>
		<b>Auto Subtract</b>	<i>0</i>
		<b>Auto SCIC</b>	<i>Off</i>
<b>CONTRAST</b>		<b>MULTI-STATION</b>	
<b>Contrast Yes/No</b>		<b>Number of Stations</b>	<i>0</i>
<i>No</i>		<b>Station Number</b>	<i>0</i>